



SINGLE-PASS FLUSHING

200EVO application

imt
a TECNIPLAST company
PHARMA

SINGLE-PASS WATER FLUSHING

Do you need to pre-wash to remove residue from your load?

Do you have low water flowrate?

Do you need a GMP single pass (non-recirculated) final rinse?

IWT can provide you with a

PHARMA GRADE BUFFER & BOOSTER WATER SYSTEM

- ▶ Variable capacity according to your needs
- ▶ Variable speed, sanitary grade pump
- ▶ Set the amount of water you need for each cycle phase
- ▶ Adjustable water pressure and temperature
- ▶ Fully drainable, with dedicated drying phase
- ▶ Self-disinfection cycle

The equipment can be used as:

- Buffer tank to shorten the cycle phase
 - Single pass pre-wash
 - Single pass GMP final rinse

Double configuration:

- ON-BOARD, to minimize machine footprint
- STAND-ALONE, on side



to find out more information:

info@iwtpharma.com | www.iwtpharma.com