

# C-LINE

High Pressure BIN Washer



**mw**<sup>®</sup>  
a TECNIPLAST company  
PHARMA



# C-LINE

## High Pressure BIN Washer

C-LINE MODEL	STANDARD HEIGHT	EXTRA HEIGHT
Chamber Footprint [mm - inch]	1400 x 1300 - 55.1 x 51.2	
Chamber Height [mm - inch]	1900 - 74.8	2350 - 92.5
Unit Footprint [mm - inch]	2025 x 1866 - 79.7 x 73.5	
Unit Height [mm - inch]	3025 - 119.1	3790 - 149.2
Sump Requirements [mm - inch]	300 - 12	

### EXCLUSIVE BENEFITS!

- ▶ Flexible and versatile: ideal to process IBCs of multiple sizes
- ▶ Made of AISI 316L (1.4404), GMP and FDA compliant
- ▶ No cross-contamination: single pass washing principle and no water recirculation
- ▶ Up to 70 Bar (1015 PSI) for great water mechanical action and quick cycles
- ▶ Reduced water and chemical consumptions: only 40 l/min (9 gal/min) water flow
- ▶ Modular design: chamber, water and drying handling units to minimize equipment footprint in classified areas
- ▶ Mobile hydrokinetic head for the cleaning of bin internal surfaces and floor-mounted hydrokinetic head for bottom valves cleaning
- ▶ Ideal for automatic high containment valves handling



to find out more information:

[info@iwtpharma.com](mailto:info@iwtpharma.com) | [www.iwtpharma.com](http://www.iwtpharma.com)

CAT MKTBCR000126P - Rev Nov 2022

IWT reserves the right to change or modify product and/or specifications without notice or obligation.