

# M-LINE

High Pressure Mobile Washer



**iwt**  
TECNIPLAST company  
PHARMA

# M-LINE


## High Pressure Mobile Washer

M-LINE	
Footprint [mm - inch]	755 x 1440 - 29.7 x 56.7
Height [mm - ft]	1250 - 49.2
Weight [kg - lbs]	300 - 660
Pressure [bar - PSI]	70 - 1015
Waterflow [l/min - gal/min]	40 - 9

### EXCLUSIVE BENEFITS!

- ▶ The ideal washing system for IBCs, drums, BIN tanks (lavabin) and, more generally, large devices in the pharmaceutical industries.
- ▶ Suitable for a wide range of applications (CIP and COP)
- ▶ Innovative design to improve ergonomics thanks to reduced weight
- ▶ Made of AISI316L, cGMP and FDA compliant
- ▶ Up to 70bar (1015 PSI) for great water mechanical action and quick cycles
- ▶ Only 40 l/min (9 gal/min) water flow for reduced water and chemical consumptions

**CHECK THE M-LINE APPLICATION POSTCARDS TO FIND OUT MORE!**



to find out more information:

[info@iwtpharma.com](mailto:info@iwtpharma.com) | [www.iwtpharma.com](http://www.iwtpharma.com)

CAT MKTBRC000120J - Rev Nov 2022

IWT reserves the right to change or modify product and/or specifications without notice or obligation.