



**NEW**





2

## COMPLETE CLEANING PROCESSES IN ONE UNIT!

CONTACT  
PARTS

IBCs



### IBCs CLEANING PROCESS:

- Telescopic cGMP Hydrokinetic header to reach any internal bin surfaces
- Bottom valve coverage via dedicated in-floor spray ball

### CONTACT PARTS CLEANING PROCESS:

- Automatic Quick-Lock Connection
- Full in-house designing and manufacturing of customized racks

### FULL cGMP DESIGN, FEATURING:

- Variable hydrodynamic control, high washing pressure and water flow tuned around every load need
- Hepa filtered drying system for bin, parts, chamber and piping

### MINIMAL BUILDING AND OPERATIONAL IMPACT:

- Only 250mm – 10" sump requirements
- Full frontal access for maintenance

### PATENT PENDING

### OZONE SANITIZATION SYSTEM

- Oxidation as a method to enhance the cleaning and sterility assurance chain
- Overall sanitization of the machine wet surfaces
- Outstanding validatable results thanks to an entirely repeatable process

