

# 2CLEAN

IBC & Contact Part Washer



**iwt**  
TECNIPLAST company  
PHARMA

# 2CLEAN

IBC & Contact Part Washer



Watch the  
2CLEAN  
in action!



## DOUBLE THE BENEFITS!

### IBCs CLEANING PROCESS

- ▶ Telescopic cGMP Hydrokinetic header to reach any internal bin surfaces
- ▶ Bottom valve coverage via dedicated in-floor spray ball

### CONTACT PARTS CLEANING PROCESS

- ▶ Automatic Quick-Lock Connection
- ▶ Full in-house designing and manufacturing of customized racks

### FULL cGMP DESIGN, FEATURING:

- ▶ Variable hydrodynamic control, high washing pressure and water flow to meet every load need
- ▶ HEPA filtered drying system for bin, parts, chamber and piping

### MINIMAL BUILDING & OPERATIONAL IMPACT:

- ▶ Only 250mm – 10” sump requirements
- ▶ Full frontal access for maintenance

2CLEAN MODEL	STANDARD	XL
Chamber Dimensions [mm - inch]	1450 x 1350 - 57.1 x 53.1	
Chamber Height [mm - inch]	2000 - 78.7	2400 - 94.4
Unit Dimensions [mm - inch]	3000 x 1600 - 118.1 x 63	
Unit total volume [m <sup>3</sup> – ft <sup>3</sup> ]	3000 - 118.1	3650 - 143.7
Ratio	250 - 10	

to find out more information:

[info@iwtpharma.com](mailto:info@iwtpharma.com) | [www.iwtpharma.com](http://www.iwtpharma.com)

CAT MKTBRC000129S - Rev Nov 2022

IWT reserves the right to change or modify product and/or specifications without notice or obligation.