





watch the 2Clean
in action!

DOUBLE THE BENEFITS!

IBCs CLEANING PROCESS:

- Telescopic cGMP Hydrokinetic header to reach any internal bin surfaces
- Bottom valve coverage via dedicated in-floor spray ball

CONTACT PARTS CLEANING PROCESS:

- Automatic Quick-Lock Connection
- Full in-house designing and manufacturing of customized racks

FULL cGMP DESIGN, FEATURING:

- Variable hydrodynamic control, high washing pressure and water flow to meet every load need
- HEPA filtered drying system for bin, parts, chamber and piping

MINIMAL BUILDING & OPERATIONAL IMPACT:

- Only 250mm – 10" sump requirements
- Full frontal access for maintenance

to find out more information: