

NEW 300 EVO

cGMP flexible Part Washer



iwt[®]
a **TECNIPLAST** company
PHARMA

NEW!

300 EVO

cGMP flexible Part Washer

Engineered to clean product contact parts from solid, semi-solid, liquid and aseptic manufacturing sites, from tableting machines to sterile filling lines

Modular design in compliance with the industry standards to allow quick turnaround projects

- ▶ Automatic sliding glass door design to fit any operational flows. Suitable for double bioseal installation
- ▶ Multiple configurations available
- ▶ Variable hydrodynamic control, high washing pressure and water flow tuned around load needs
- ▶ Up to 30 fully customizable washing programs from pre-washing to drying and cooling
- ▶ PanelPC configuration for remote connectivity and **INDUSTRY 4.0 integration**; cutting edge automation system GAMP5 to guarantee full data traceability and storage
- ▶ **Product contact parts entirely made of AISI316L (1.4404)**, roughness < 0.5 µm ASME-BPE compliant
- ▶ Multiple water inlet and chemical dosing systems; single Pass Rinse feature via dedicated tank; HEPA filtered drying system
- ▶ **Full in-house design and manufacturing of customized racks** connected via automatic quick-lock patented system
- ▶ **Comprehensive set of handling solutions** for manoeuvring, stocking and loading/unloading
- ▶ **Hinged process pump for ease of maintenance**

300 EVO MODEL	100	140	200
Chamber Nominal Volume [m ³ – ft ³]	> 1.0 – 35.3	> 1.4 – 49.4	> 2.0 – 70.6
Footprint [m ² – ft ²]	~ 2.8 – 30	~ 3.6 – 38.8	~ 4.9 – 52.7
Unit height [m – ft]	2.7 – 8.9		

Three different sizes (1000, 1400, 2000 liters - 265, 370, 530 gal) to maximize and optimize cleaning batches



to find out more information:

info@iwtpharma.com | www.iwtpharma.com

CAT MKTBRC000131M - Rev Nov 2022

IWT reserves the right to change or modify product and/or specifications without notice or obligation.