


**BLENDER
CLEANING**

BLENDER CLEANING

Do you have to...

- Automatically clean a large and tall blender
- Reduce the down-time of the equipment increasing productivity
- Implement a fully validated cleaning process, faster and cost effective
 - Increase the safety for the operators



IWT have the perfect cleaning solution: **M-LINE HIGH PRESSURE WASHER**

- Multiple hydrokinetic lances connection for total internal volume coverage
- 70Bar | 1015Psi water pressure for an impressive scraping effect on surfaces
 - 7 min total cycle time
 - Low water consumption: only 40l/min | 10gal/min
- - 82% water, - 90% time compared to previous manual cleaning process

to find out more information:

info@iwtpharma.com | www.iwtpharma.com