




C-LINE



watch the C-Line
in action!

EXCLUSIVE BENEFITS!

- Flexible and versatile: ideal to process IBCs of multiple sizes
 - Made of AISI 316L (1.4404), GMP and FDA compliant
 - No cross-contamination: single pass washing principle and no water recirculation
- Up to 70 Bar (1015 PSI) for great water mechanical action and quick cycles
 - Reduced water and chemical consumptions: only 40 l/min (9 gal/min) water flow
- Modular design: chamber, water and drying handling units to minimize equipment footprint in classified areas
- Mobile hydrokinetic head for the cleaning of bin internal surfaces and floor-mounted hydrokinetic head for bottom valves cleaning
- Ideal for automatic high containment valves handling

to find out more information: