



FILLING

LIQUID FORMS

200EVO application

iwat[®]
TECNIPLAST company
PHARMA

FILLING

LIQUID FORMS

Do you have to clean aseptic filling line parts as:

Dosing pumps | Filling needles |
Connection hoses and fittings |
Stopper Bowls |
Hoppers and other accessories



IWT have the perfect cleaning solution:

200EVO cGMP CONTACT PARTS WASHER

- ▶ Customized loading baskets with “punctual” cleaning to ensure perfect coverage of the load
- ▶ Dedicated spraying orifices to flush filling needles and connection hoses/fittings
- ▶ Direct WFI final rinse to provide the highest level of hygiene
- ▶ Final rinse in-line conductivity and TOC monitoring for high cleaning quality assurance
- ▶ Perfect drying facility with HEPA filtered air to ensure no water residues thus preventing any microbial growth
- ▶ Hard and electronic batch report (21CFR Part11 compliant) to guarantee full data traceability and storage

to find out more information:

info@iwtpharma.com | www.iwtpharma.com