




2CLEAN
cGMP IBC & CONTACT PART WASHER

NEW



CLEANING



2

**COMPLETE
CLEANING PROCESSES
IN ONE UNIT!**

CONTACT
PARTS

IBCs



IBC's CLEANING PROCESS:

- Telescopic GMP Hydrokinetic header to reach any internal bin surfaces
- Bottom valve coverage via dedicated in-floor spray ball

CONTACT PARTS CLEANING PROCESS:

- Automatic Quick-Lock Connection
- Full in-house designing and manufacturing of customized racks

FULL CGMP DESIGN, FEATURING:

- Variable hydrodynamic control, high washing pressure and water flow tuned around every load need
- Hepa filtered drying system for bin, parts, chamber and piping


MINIMAL BUILDING AND OPERATIONAL IMPACT:

- Only 250mm – 10" sump requirements
- Full frontal access for maintenance

PATENT PENDING

OZONE SANITIZATION SYSTEM

- Oxidation as a method to enhance the cleaning and sterility assurance chain
- Overall sanitization of the machine wet surfaces
- Outstanding validatable results thanks to an entirely repeatable process



COMPANY WITH
QUALITY SYSTEM
CERTIFIED BY DNV GL
= ISO 9001 =

